
UV-1800+UVProbe

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Photometric : (Multi λ 多段波長分析)

(1) 按下 進入 **Photometric Method Wizard** (設定所需波長)

Photometric Method Wizard - [Wavelengths]

Wavelength Type: Point Column Name: WL800.0

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Photometric Method Wizard - [Calibration]

Type: Multi Point Column Name:

Formula: Fixed Wavelength Units: mg/l

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
Sample Table - [Active]					
	Sample ID	Type	Ex	WL800.0	Comments
1 *	TEST1				
2	TEST2				
3					

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
Calibration 畫面

Photometric Method Wizard - [Calibration]

Type: **Multi Point** Column Name:

Formula: **Fixed Wavelength** Units: **mg/l**

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


Photometric Method Wizard - [Measurement Parameters(Standard)]

STD

Data Acquired By: User Entry Instrument

Sample Repetitions: Prompt before repeat

Delay Sample Read Time: seconds

Auto Scan Timed Repeat None

Data Acquired By: User Entry Instrument :

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Sample Repetitions: Prompt before repeat :

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

The image shows a dialog box titled "Photometric Method Wizard - [Measurement Parameters(Sample)]". It features an "SMP" logo in the top left. The "Data Acquired By:" section has two radio buttons: "User Entry" (unselected) and "Instrument" (selected). The "Sample Repetitions:" section has a text box containing "1" and a "Prompt before repeat" checkbox (unchecked). The "Delay Sample Read" checkbox is unchecked, and the "Time:" section has a text box containing "0" and a dropdown menu set to "seconds". At the bottom, there are three radio buttons: "Auto Scan" (unselected), "Timed Repeat" (unselected), and "None" (selected).

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The image shows a dialog box titled "Photometric Method Wizard - [File Properties]". It contains an information icon and a message: "When you create a Method file or a Standard file (not a Data file), please click the Finish button now! You can edit some properties of this page later in File > Properties > Summary Tab." Below this, there are fields for "Filename:" (containing "C:\Uvpc39\DATA\File_040531_161747.pho"), "Title:", "Analyst:", and "Comments:". At the bottom, there is a note: "Method Property Sheet will be displayed when the Finish button is clicked. If you need to use advanced parameters of Photometric such as Equations, Pass/Fail, Instrument and Attachments. Please set up these pages." and three buttons: "< 上一步(B)", "完成", and "取消".

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Sample Table								
	Sample ID	Type	Ex	Conc	WL580.00	WL550.00	Result	Comments
14	14	Unknown		24.077	0.537	0.301	0.236	imported from old data file
15	15	Unknown		22.813	0.479	0.255	0.224	imported from old data file
16	16	Unknown		7.190	0.585	0.513	0.072	imported from old data file
17	17	Unknown		0.779	0.163	0.153	0.010	imported from old data file
18	18	Unknown		6.777	0.346	0.278	0.068	imported from old data file
19	19	Unknown		30.866	0.654	0.352	0.302	imported from old data file
20	20	Unknown		48.601	0.795	0.327	0.468	imported from old data file
21	21	Unknown		69.149	0.935	0.289	0.646	imported from old data file
22	22	Unknown		3.741	0.147	0.109	0.039	imported from old data file

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Step 2 – 連結 Peak Pick Table 到 Active Spectrum

The screenshot displays the IVProba software interface. At the top, there is a menu bar with options: File, Edit, View, Insert, Tools, Window, Help. Below the menu bar is a toolbar with various icons. The main workspace is divided into several sections:

- Spectrum Plot:** Located in the upper left, it shows a plot with a y-axis ranging from 0.00 to 0.20 and an x-axis labeled 'nm.' ranging from 0.00 to 1.00. The plot area is currently empty.
- Peak Pick Table:** A table with the following columns: No., P/V, Wavelength nm., ..., and Description. The table is currently empty.
- Peak Pick Properties Dialog:** A dialog box titled 'Peak Pick Properties' is open in the lower right. It has tabs for 'Data', 'Pinning Options', 'General', 'Peaks', and 'Values'. The 'Data' tab is selected, showing a list of spectra with checkboxes:
 - Active Spectrum
 - Most Recent Spectrum 1
 - Most Recent Spectrum 2
 - Most Recent Spectrum 3
 - Most Recent Spectrum 4
 - Most Recent Spectrum 5
 - Most Recent Spectrum 6
 - Most Recent Spectrum 7A 'View Links' button is located at the bottom right of the dialog.

The status bar at the bottom of the window shows 'Page 1 / 1', 'Pos: 5.85, 7.13 n', and 'NUM'.

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